

ANTIMICROBIAL
Treated to Help Keep Swabs Fresher

Q-tips® Cotton Swabs

More soft
cotton at the tip* than
any other swab



Q-tips®

WARNING: Do not insert swab into ear canal. Entering the ear canal could cause injury. If used to clean ears, stroke swab gently around the

outer surface of the ear only. Keep out of reach of children. Q-tips® Cotton Swabs are biodegradable and packaging is recyclable†.

Questions: Call 1-800-265-7964
Consumer Info Center

For more helpful hints
visit www.qtips.com
*From the end of the stick
to the top of the swab.



83136288

This seal signifies that Q-tips® Cotton Swabs are made with 100% pure cotton. The seal of cotton is a registered service mark/trademark of Cotton Incorporated.



© UNILEVER
Trumbull, CT 06611
Made in the U.S.A.

†Packaging is coded for recycling. Check to see if facilities exist in your area.

300 COTTON SWABS

ANTIMICROBIAL
Treated to Help Keep Swabs Fresher
Q-tips®
Cotton Swabs

300 COTTON SWABS

ANTIMICROBIAL
Treated to Help Keep Swabs Fresher

Q-tips®
Cotton Swabs

300 COTTON SWABS